



Product Datasheet

PKCδ (Phospho-Ser645) Antibody (orb14972)

Description	PKCδ (Phospho-Ser645) Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q05655
Expiration Date	12 months from date of receipt.



Western blot analysis of PKCδ(Phospho-Ser645) in PC12 cells treated with PMA.



Immunohistochemical analysis of formalin-fixed tissue showing PKCδ (Phospho-Ser645) staining.