


Product Datasheet

Smad2 (Phospho-Thr220) Antibody (orb14987)

Description	Smad2 (Phospho-Thr220) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q15796
Expiration Date	12 months from date of receipt.



Western blot analysis of Smad2(Phospho-T...