


Product Datasheet

Tyrosine Hydroxylase (Phospho-Ser19) Antibody (orb14992)

Description	Tyrosine Hydroxylase (Phospho-Ser19) Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P07101
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...