
Product Datasheet

Biotinylated Human GLP1R Protein, His, (orb1499256)

Description

Biotinylated Human GLP1R, His, is

(HEK293). It contains AA Ala 21 - G...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

GLP1R

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 17.3 kDa. The protein migrates as 30-34 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Biotinylated Human GLP1R Protein, His, (orb1499256) is expressed from human 293 cells (HEK293). It contains AA Ala 21 - Glu 139 (Accession # P43220-1).

MW

17.3 kDa

Uniprot ID[P43220](#)**Expiration Date**

6 months from date of receipt.