


Product Datasheet

FITC-Labeled Human IL-1 RAcP / IL-1 R3 Protein, His Tag (orb1499259)

Description	FITC-Labeled Human IL-1 RAcP Protein, His Tag is expressed from human 293 cells (HEK293). It contain...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	FITC
Target	IL-1 RAcP
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-8×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.69 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	FITC-Labeled Human IL-1 RAcP Protein, His Tag (orb1499259) is expressed from human 293 cells (HEK293). It contains AA 21 Ser - 359 Glu (Accession # Q9NPH3-1).
MW	39.69 kDa
Expiration Date	6 months from date of receipt.



Immobilized
Human
Anti-IL1RAP
Human
IgG1...