

Product Datasheet

IKB Alpha (phospho-Ser32/Ser36) antibody (orb1499271)

Description	IKB Alpha (phospho-Ser32/Ser36) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human phospho-IKB Alpha (Ser32 / Ser36)
Target	phospho-IKB Alpha (Ser32 / Ser36)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	35 kDa
Uniprot ID	P25963
Expiration Date	12 months from date of receipt.