

Product Datasheet

eIF4EBP1 (phospho-Thr37/Thr46) antibody (orb1499304)

Description	eIF4EBP1 (phospho-Thr37/Thr46) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human phospho-eIF4EBP1 (Thr37 / Thr46)
Target	phospho-eIF4EBP1 (Thr37 / Thr46)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	13 kDa
Uniprot ID	Q13541
Expiration Date	12 months from date of receipt.