

Product Datasheet

SIRT1 (phospho-Ser27) antibody (orb1499327)

Description	SIRT1 (phospho-Ser27) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human phospho-SIRT1 (Ser27)
Target	phospho-SIRT1 (Ser27)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	81 kDa
Uniprot ID	Q96EB6
Expiration Date	12 months from date of receipt.