

Product Datasheet

HCMV pp65 antibody (orb1499414)

Description	HCMV pp65 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant HCMV pp65 protein
Target	HCMV pp65
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	65 kDa
Uniprot ID	Q6SW59
Expiration Date	12 months from date of receipt.