



Product Datasheet

Her3/ErbB3 (phospho-Tyr1328) Antibody (orb14999)


Description	Her3/ErbB3 (phospho-Tyr1328) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P21860
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of M...



Immunohistochemical analysis of formalin...



Western blot analysis of Her3/ErbB3(phos...