



Product Datasheet

HDAC4/HDAC5/HDAC9 (phospho-Ser246/259/220) Antibody (orb15003)

Description	HDAC4/HDAC5/HDAC9 (phospho-Ser246/259/220) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P56524
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of HDAC4/HDAC5/HDA...