


Product Datasheet


MARCKS (phospho-Ser170) Antibody (orb15008)

Description	MARCKS (phospho-Ser170) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P29966
Expiration Date	12 months from date of receipt.



Western blot analysis of MARCKS(phospho-

...



Immunohistochemical analysis of formalin...