

## Product Datasheet

# DPYSL2 (phospho-Thr514+Ser518) Antibody Blocking Peptide (orb1501342)

<b>Description</b>	DPYSL2 (phospho-Thr514+Ser518) Antibody Blocking Peptide
<b>Target</b>	DPYSL2
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q16555</a>
<b>Expiration Date</b>	12 months from date of receipt.