

Product Datasheet

Eph receptor A2+A3+A4 (phospho-Tyr588 + Tyr596) Antibody Blocking Peptide (orb1502007)

Description	Eph receptor A2+A3+A4 (phospho-Tyr588 + Tyr596) Antibody Blocking Peptide
Target	EPHA2
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Uniprot ID	P29317
Expiration Date	12 months from date of receipt.