

## Product Datasheet

### STK26 Antibody Blocking Peptide (orb1502512)

|                        |  |
|------------------------|--|
| <b>Description</b>     | STK26 Antibody Blocking Peptide  |
| <b>Target</b>          | STK26  |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <a href="#">Q9P289</a>   |
| <b>Expiration Date</b> | 12 months from date of receipt.  |