


Product Datasheet

EGFR (phospho-Tyr1197) Goat Polyclonal Antibody (orb15029)

Description	EGFR (phospho-Tyr1197) Goat Polyclonal Antibody
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.



Western blot analysis of EGFR(phospho-Ty...