

## Product Datasheet

# SCNN1B (phospho-Ser633) Antibody Blocking Peptide (orb1503409)

<b>Description</b>	SCNN1B (phospho-Ser633) Antibody Blocking Peptide
<b>Target</b>	SCNN1B
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P51168</b>
<b>Expiration Date</b>	12 months from date of receipt.