

## Product Datasheet

# P53 (phospho-Thr55) Antibody Blocking Peptide (orb1503479)

|                        |  |
|------------------------|--|
| <b>Description</b>     | P53 (phospho-Thr55) Antibody<br>Blocking Peptide                                       |
| <b>Target</b>          | TP53   |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for<br>one year. Avoid repeated<br>freeze/thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>P04637</b>  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |