

## Product Datasheet

# ELK1 (phospho-Ser389) Antibody Blocking Peptide (orb1503688)

|                        |  |
|------------------------|--|
| <b>Description</b>     | ELK1 (phospho-Ser389) Antibody<br>Blocking Peptide                                     |
| <b>Target</b>          | ELK1   |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for<br>one year. Avoid repeated<br>freeze/thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>P19419</b>  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |