

Product Datasheet

AKT1+AKT2 (phospho-Thr308+Thr309) Antibody Blocking Peptide (orb1503798)

Description	AKT1+AKT2 (phospho-Thr308+Thr309) Antibody Blocking Peptide
Target	AKT1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Uniprot ID	P31749
Expiration Date	12 months from date of receipt.