

## Product Datasheet

# H5N1 Nucleoprotein Antibody Blocking Peptide (orb1503971)

|                        |                                                                                        |
|------------------------|----------------------------------------------------------------------------------------|
| <b>Description</b>     | H5N1 Nucleoprotein Antibody<br>Blocking Peptide                                        |
| <b>Target</b>          | H1N1 N                                                                                 |
| <b>Form/Appearance</b> | Lyophilized powder                                                                     |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for<br>one year. Avoid repeated<br>freeze/thaw cycles. |
| <b>Note</b>            | For research use only                                                                  |
| <b>Uniprot ID</b>      | <b>B4URE0</b>                                                                          |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                        |