

Product Datasheet

GSK3 Alpha/Beta (phospho-Tyr279/Tyr216) Antibody Blocking Peptide (orb1504533)

| | |
|------------------------|--|
| Description | GSK3 Alpha/Beta (phospho-Tyr279/Tyr216) Antibody Blocking Peptide |
| Target | GSK3B |
| Storage | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Uniprot ID | P49841 |
| Expiration Date | 12 months from date of receipt. |