

## Product Datasheet

# Stat4 (phospho-Tyr693) Antibody Blocking Peptide (orb1505089)

<b>Description</b>	Stat4 (phospho-Tyr693) Antibody Blocking Peptide
<b>Target</b>	STAT4
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	KG(p-Y)VP
<b>Uniprot ID</b>	<a href="#">Q14765</a>
<b>Expiration Date</b>	12 months from date of receipt.