

## Product Datasheet

# SGK1 (phospho-Ser78) Antibody Blocking Peptide (orb1505123)

|                        |  |
|------------------------|--|
| <b>Description</b>     | SGK1 (phospho-Ser78) Antibody<br>Blocking Peptide                                      |
| <b>Target</b>          | SGK1   |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for<br>one year. Avoid repeated<br>freeze/thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <a href="#">O00141</a>   |
| <b>Expiration Date</b> | 12 months from date of receipt.  |