

## Product Datasheet

# PDGFRB (phospho-Tyr1009) Antibody Blocking Peptide (orb1505231)

<b>Description</b>	PDGFRB (phospho-Tyr1009) Antibody Blocking Peptide
<b>Target</b>	PDGFRB
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P09619</b>
<b>Expiration Date</b>	12 months from date of receipt.