

## Product Datasheet

# PDGFRA (phospho-Tyr849) /PDGF Receptor beta (phospho-Tyr857) Antibody Blocking Peptide (orb1505232)

<b>Description</b>	PDGFRA (phospho-Tyr849) /PDGF Receptor beta (phospho-Tyr857) Antibody Blocking Peptide
<b>Target</b>	PDGFRA
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P16234</b>
<b>Expiration Date</b>	12 months from date of receipt.