

Product Datasheet

PAK1 (phospho-Ser199/204) /PAK2 (phospho-Ser192/197) Antibody Blocking Peptide (orb1505243)

Description	PAK1 (phospho-Ser199/204) /PAK2 (phospho-Ser192/197) Antibody Blocking Peptide
Target	PAK1
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Uniprot ID	Q13153
Expiration Date	12 months from date of receipt.