

Product Datasheet

STK4 (phospho-Thr183) /Mst2 (phospho-Thr180) ; STK4/3 (phospho-Thr183/180) Antibody Blocking Peptide (orb1505294)

Description	STK4 (phospho-Thr183) /Mst2 (phospho-Thr180) ; STK4/3 (phospho-Thr183/180) Antibody Blocking Peptide
Target	STK4
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Uniprot ID	Q13043
Expiration Date	12 months from date of receipt.