

Product Datasheet

PPP1R1B (phospho-Thr34) Antibody Blocking Peptide (orb1505435)

| | |
|------------------------|--|
| Description | PPP1R1B (phospho-Thr34) Antibody Blocking Peptide |
| Target | PPP1R1B |
| Storage | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Uniprot ID | Q9UD71 |
| Expiration Date | 12 months from date of receipt. |