

Product Datasheet

ATP1A1 (phospho-Ser23) Antibody Blocking Peptide (orb1505478)

| | |
|------------------------|----------------------------------------------------------------------------------|
| Description | ATP1A1 (phospho-Ser23) Antibody Blocking Peptide |
| Target | ATP1A1 |
| Storage | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Uniprot ID | P05023 |
| Expiration Date | 12 months from date of receipt. |