

## Product Datasheet

# ACLY (phospho-Ser455) Antibody Blocking Peptide (orb1505494)

|                        |                                                                                  |
|------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>     | ACLY (phospho-Ser455) Antibody Blocking Peptide                                  |
| <b>Target</b>          | ACLY                                                                             |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only                                                            |
| <b>Uniprot ID</b>      | <b>P53396</b>                                                                    |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                  |