

Product Datasheet

ATP1A1 (phospho-Ser16) ; Na,K-ATPase alpha-1 (phospho-Phospho Ser16) Antibody Blocking Peptide (orb1505496)

Description	ATP1A1 (phospho-Ser16) ; Na,K-ATPase alpha-1 (phospho-Phospho Ser16) Antibody Blocking Peptide
Target	ATP1A1
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Uniprot ID	P05023
Expiration Date	12 months from date of receipt.