

Product Datasheet

p44/42 MAPK (phospho-ERK1/2) (phospho-Thr202/Tyr204) Antibody Blocking Peptide (orb1505512)

Description	p44/42 MAPK (phospho-ERK1/2) (phospho-Thr202/Tyr204) Antibody Blocking Peptide
Target	MAPK3
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Uniprot ID	Q63844
Expiration Date	12 months from date of receipt.