

Product Datasheet

beta 1 Adrenergic Receptor antibody (FITC) (orb15064)

Description	beta 1 Adrenergic Receptor antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide (please contact us for details)
Target	beta 1 Adrenergic Receptor (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
Uniprot ID	P18090, P08588
Dilution Range	IF/ICC:1:50-200 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.