
Product Datasheet

FGFA/FHFA (pan) antibody (orb150651)

Description

Mouse monoclonal to FGFA (Biotin).

called FHF2 is a protein...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Biotin

Tested

ICC, IHC, WB

Applications
Immunogen

Synthetic peptide amino acids 2-18 (AAAIASSLIRQKRQARE) of human FHF2A. 100% identical to rat, 94% identical to mouse. > 80% identity with FGF12A/FHF1A, FGF14A/FHF4A and FGF11A/FHF3A.

Target

FGF13

Preservatives

PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-448 was sufficient for detection of FGFA/FHFA (pan) in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2b

Clonality

Monoclonal

MW

30kDa

Uniprot ID

[Q92913](#)

NCBI

[AAH34340](#)

Entrez

[2258](#)

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.