

---

## Product Datasheet

### FGFA/FHFA (pan) antibody (orb150655)

**Description**

Mouse monoclonal to FGFA (PerCP).

called FHF2 is a protein ...

**Species/Host**

Mouse

**Reactivity**

Human, Mouse, Rat

**Conjugation**

PerCP

**Tested Applications**

ICC, IHC, WB

**Immunogen**

Synthetic peptide amino acids 2-18 (AAAIASSLIRQKRQARE) of human FHF2A. 100% identical to rat, 94% identical to mouse. > 80% identity with FGF12A/FHF1A, FGF14A/FHF4A and FGF11A/FHF3A.

**Target**

FGF13

**Preservatives**

PBS pH 7.4, 50% glycerol, 0.1% sodium azide \*Storage buffer may change when conjugated

**Concentration**

1 mg/ml

**Storage**

Conjugated antibodies should be stored at 4°C

**Note**

For research use only

**Application notes**

1 µg/ml of SMC-448 was sufficient for detection of FGFA/FHFA (pan) in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

**Isotype**

IgG2b

**Clonality**

Monoclonal

**MW**

30kDa

**Uniprot ID**

[Q92913](#)

**NCBI**

[AAH34340](#)

**Entrez**

[2258](#)

**Dilution Range**

WB (1:1000)

**Expiration Date**

12 months from date of receipt.