

Product Datasheet

HCMV UL49 (Human Cytomegalovirus UL49 Gene; Human herpesvirus 5.□ Antibody Blocking Peptide (orb1506609)

Description	HCMV UL49 (Human Cytomegalovirus UL49 Gene; Human herpesvirus 5.□ Antibody Blocking Peptide
Target	UL49
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Uniprot ID	P16786
Expiration Date	12 months from date of receipt.