

## Product Datasheet

### T0901317 (orb15069)

<b>Description</b>	T0901317 is a potent, high affinity liver X receptor (LXR) agonist (EC50 ~ 50 nM, Kd values are 7 an...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{17}H_{12}F_9NO_3S$
<b>MW</b>	481.33
<b>CAS Number</b>	[293754-55-9]
<b>SMILES</b>	<chem>C(S(N(CC(F)(F)F)C1=CC=C(C(O)(C(F)(F)F)C(F)(F)F)C=C1)(=O)=O)1=CC=CC=C1</chem>
<b>Expiration Date</b>	12 months from date of receipt.