

## Product Datasheet

# BCL2L1 (phospho-Ser62) Antibody Blocking Peptide (orb1507584)

<b>Description</b>	BCL2L1 (phospho-Ser62) Antibody Blocking Peptide
<b>Target</b>	BCL2L1
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q07817</a>
<b>Expiration Date</b>	12 months from date of receipt.