

## Product Datasheet

### NLRP14 Antibody Blocking Peptide (orb1508447)

|                        |  |
|------------------------|--|
| <b>Description</b>     | NLRP14 Antibody Blocking Peptide   |
| <b>Target</b>          | NLRP14   |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <a href="#">Q86W24</a>   |
| <b>Expiration Date</b> | 12 months from date of receipt.  |