

## Product Datasheet

### TRANK1 Antibody Blocking Peptide (orb1510898)

|                        |                                                                                  |
|------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>     | TRANK1 Antibody Blocking Peptide                                                 |
| <b>Target</b>          | TRANK1                                                                           |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only                                                            |
| <b>Uniprot ID</b>      | <a href="#">O15050</a>                                                           |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                  |