

## Product Datasheet

### CHTOP Antibody Blocking Peptide (orb1511945)

|                        |                                                                                  |
|------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>     | CHTOP Antibody Blocking Peptide                                                  |
| <b>Target</b>          | CHTOP                                                                            |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only                                                            |
| <b>Uniprot ID</b>      | <b>Q9Y3Y2</b>                                                                    |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                  |