

## Product Datasheet

# 14-3-3 Tau (phospho-Ser232) Antibody Blocking Peptide (orb1512969)

<b>Description</b>	14-3-3 Tau (phospho-Ser232) Antibody Blocking Peptide
<b>Target</b>	YWHAQ
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P27348</b>
<b>Expiration Date</b>	12 months from date of receipt.