

## Product Datasheet

# NR3C1 (phospho-Ser226) Antibody Blocking Peptide (orb1513294)

|                        |  |
|------------------------|--|
| <b>Description</b>     | NR3C1 (phospho-Ser226) Antibody Blocking Peptide                                 |
| <b>Target</b>          | NR3C1  |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>P04150</b>  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |