

## Product Datasheet

# GATA2 (phospho-Ser401) Antibody Blocking Peptide (orb1513377)

<b>Description</b>	GATA2 (phospho-Ser401) Antibody Blocking Peptide
<b>Target</b>	GATA2
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P23769</b>
<b>Expiration Date</b>	12 months from date of receipt.