

## Product Datasheet

# PLB (phospho-phospho Ser16+Thr17) Antibody Blocking Peptide (orb1515462)

|                        |  |
|------------------------|--|
| <b>Description</b>     | PLB (phospho-phospho Ser16+Thr17) Antibody Blocking Peptide                      |
| <b>Target</b>          | PLN  |
| <b>Storage</b>         | Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>P26678</b>  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |