

Product Datasheet

MST4 + MST3 + STK25 (phospho-Thr178 + Thr190 + Thr174) Antibody Blocking Peptide (orb1515784)

Description	MST4 + MST3 + STK25 (phospho-Thr178 + Thr190 + Thr174) Antibody Blocking Peptide
Target	STK25
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Uniprot ID	O00506
Expiration Date	12 months from date of receipt.