

Product Datasheet

Human CD19 Protein (FITC) (orb1516290)

Description	Recombinant Human CD19 Protein conjugated with FITC
Endotoxins	< 1.0 EU/μg as determined by LAL
Target	CD19
Preservatives	Supplied as 0.22μm filtered solution in PBS (pH=7.4)
Form/Appearance	Liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	20-291/556
Purity	> 95% as determined by Tris-Bis PAGE
Uniprot ID	P15391
Expiration Date	6 months from date of receipt.