

## Product Datasheet

### Human DR4 Protein (orb1516308)

<b>Description</b>	Recombinant Human DR4 Protein
<b>Endotoxins</b>	< 1.0 EU/μg as determined by LAL
<b>Target</b>	TNFRSF10A
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	C-His-Avi
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	34-239/468
<b>Purity</b>	> 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC
<b>MW</b>	Due to glycosylation, the protein migrates to 25-30 kDa based on Tris-Bis PAGE result. kDa
<b>Uniprot ID</b>	<a href="#">O00220</a>
<b>Expiration Date</b>	6 months from date of receipt.