

Product Datasheet

Human Mature TGF Beta 2 Protein (orb1516320)

Description	Recombinant Human Mature TGF Beta 2 Protein
Endotoxins	< 1.0 EU/μg as determined by LAL
Target	TGFB2
Form/Appearance	Lyophilized powder
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in 4 mM HCl.
Protein Sequence	303-414/414
Purity	> 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC
MW	Due to glycosylation, the protein migrates to 13-15 kDa based on Tris-Bis PAGE result. kDa
Uniprot ID	P61812
Expiration Date	6 months from date of receipt.